VP 500-L ...

**Vötsch**  
Industrietechnik

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# ... Photostability Test Chamber VP 500-L

## The Application ...

One of the most important requirements with regard to photostability testing is homogeneous illumination of the specimens.

For this reason and on account of the fact that the illumination of the light decreases with the distance, it is indispensable that all specimens are stored at the same distance to the illumination source. Therefore, the distance to the source of the light is constant throughout the entire storage area.

In addition to that the inhomogeneous emission of light by fluorescent tubes is balanced by special light and UV filter systems, hence homogeneous illumination of the entire storage area and, therefore, also of all specimens is achieved.

In order to fulfil the requirements by GMP or by authorities, such as e.g. the FDA, comprehensive documents are available for qualification incl. FAT, DQ, IQ, OQ and parts of PQ upon request.

## The Test Chamber...

The Vötsch **photostability test chamber VP 500-L** has an ideal light, temperature and humidity distribution and is thus able to simulate reproducible light and climatic conditions.

The lighting equipment complies with the **ICH Guideline Q1B Option 2** and enables photostability tests to be performed in less than 100 hours. All specimens are positioned at the same distance from the light source and **are evenly irradiated by a special light/UV filter system.**

This system can be equipped with appropriate light and UV sensors for recording irradiation values, e.g. the total number of lux hours as well as the total irradiation.

The required relative air humidity is generated with a patented steam humidification system (Sterile Steam System) in which the water introduced into the system is evaporated at 140°C. The humidity measurement is carried out by means of a capacitive humidity sensor.

Calibrations and spatial measurements on site as well as processing of the qualification documents and DKD (German Calibration Service) calibrations at the factory and on site are available.

Documentation of the measured values is carried out with the help of our **SIMPATI\* Pharma** software package which facilitates documentation corresponding to **FDA 21 CFR Part 11** either with the help of the control loop sensors or with the help of independent sensors.

Moreover, different analog and digital recorders (also in conformity with FDA 21 CFR Part 11) are available.

## Standard equipment

- Microprocessor monitoring and control unit SIMCON/32\*-NET with terminal MINCONTROL
- 2 shelves illuminated with UV-light and 2 shelves illuminated with white light
- Timer for light and UV-light
- Light and UV filters for ideal distribution
- Independent adjustable temperature limiter  $t_{\min} / t_{\max}$
- Interface RS 232
- Water tank 19 l, with manual and automatic water supply possibility
- Patented vapour humidification (Sterile Steam System)
- Air-cooled refrigeration unit
- Capacitive humidity measuring sensor
- Lockable door
- Door switch
- Mobile design on wheels
- Calibration of 2 temperature and 2 humidity values

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*Performance values refer to +25 °C ambient temperature. <sup>1)</sup> In accordance with IEC 60068-3-5.*

*<sup>2)</sup> Free field, 1 m distance from the front, as per DIN 45635, part 1, accuracy class 2.*

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## Important Options

- Colour touchpanel
- UV and Lux-sensors with integration function
- Analog line recorder
- Digital line recorder
- TCP/IP interface
- Networking (RS 485 interface)
- Software SIMPATI\*-**Pharma**
- Qualification documentation
- DKD calibrations
- Mapping of light- and UV-distribution
- Spatial measurements for temperature and humidity
- Independent temperature and humidity sensors for documentation purposes
- Water-cooled refrigeration unit
- Maintenance contracts

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**We reserve the right of changes in construction resulting from technical progress. Some of the illustrated systems contain optional extras.**

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**Vötsch Industrietechnik GmbH**